CLI will run, but other programs such netconfig, datalog, flash etc. may experience an error: MSVCR110.dll missing or MSVCP140.dll missing

ILTs application requires Visual C++ 2012 redistributables. 2010 and 2015 will not work with Datalight III

MSVC ####dll is the Microsoft Visual C++ Redistributable dll that is needed for projects built with

Visual Studio 2011. The **dll** letters spell this out.

MS = Microsoft, V = Visual, C = C++, R = Redistributable, P Redistributable Package.

This error appears when you wish to run a software which require the Microsoft Visual C++ Redistributable 2012. The redistributable can easily be downloaded on the Microsoft website as x86 or x64 edition

x64 edition.

Fix:

If re-installation does not fix the problem it is most likely the Visual C++ redistributables (CLI will run without the redistributables).

Older versions of Datalight software, ran with Visual C++ 2010 redistributables. Datalight III needs Visual C++ **2012** redistributables.

The install kit should fix this automatically, however, sometimes a configuration does not allow it.

Download and install the redistributables here: <u>http://www.microsoft.com/en-us/download/details.aspx?id=30679</u>

There will be a choice to install either x64 or x86, ILT recommends installing both (the computer may be 64-bit system, the ILT Apps run as x86).